



Model
1024

Evolution Series Rehab Shower/Commode Chair Model # 1024 Adjustable Seat Angle

This model introduces an easily adjustable seat slope by merely changing the rear axle location. Simply locate the Quick Release Axles and place in any of several positions to obtain a slope up to 2" front to rear. No tools are required. Integral user positioning is addressed through combining the seat angle adjustment with a multiple choice of padded seats, contoured or flat. The slotted seat opening can be moved to any direction easily and without tools. Choose from the many additional options available to achieve a tailor-made fit. All these easy adjustments make it comfortable for the user and easy for the caregiver.



Model 1024

EVOLUTION CHAIR ACCESSORIES

- | | | |
|--------------|-------|---|
| Option Codes | • A1 | Arm Troughs |
| | • A4 | *Height Adjustable Flip Back Padded Arms |
| | • C7 | *Comfortuff™ Contoured Seat (Front Open Only) |
| | • C8 | Ensolite U-Shaped Seat 2" Waterfall (4-Way Direction) |
| | • D9 | Swing Away Adjustable Laterals |
| | • D8 | *Sling Back with Stroller Handles |
| | • D10 | *Solid Adjustable Depth Back |
| | • D11 | *Slack Adjustable Depth Back |
| | • F | Seat Belt (hook and loop) |
| | • F2 | Bodypoint® Trunk Belt |
| | • G9 | *Height Adjustable Swing Away Footrests – Removable |
| | • H2 | PVC Dipped Handrings |
| | • H5 | Pneumatic Tires |
| | • H7 | Solid Treaded 24" Tire |
| | • M | Heel Loops |
| | • N | Wheel Lock Handle Extensions |
| | • P | Pan Hanger and Pan |

(*See Price List order form for option photos)

“The Gold Standard” since 1965

ACTIVEAID, Inc. 101 Activeaid Road, Box 359 Redwood Falls, Minnesota 56283-0359, U.S.A.
Call: 507-644-2951 Toll Free: 1-800-533-5330 FAX: 507-644-2468 www.activeaid.com

EVOLUTION SERIES



**Evolution Series Rehab Shower/Commode Chair
Model # 1024
Adjustable Seat Angle**

Seat width	18"
Seat depth	18"
Floor to top of seat height (front of seat)	23"
Back height (from seat)	19"
Arm height (from seat – to top of standard arm pad)	8"
Under seat clearance	18 ⁵ / ₈ ", 19 ⁵ / ₈ , 20 ⁵ / ₈ "
Distance between Quick Release Axles	15"
Footrest height adjustment	15 1/2" – 19 1/2"
Overall width	25 1/4"
Seat width between arms	20 1/2" std. 18 1/2" adj.
Overall Length	34 3/4"
Overall Length with G9, G10 (Ht. Adjustable Footrests)	43 1/2"



1024 shown with option G9
And set at lowest angle



1024 shown with options
A4, C7, D8, F2, G9, & P

Visit our Website: www.activeaid.com

“The Gold Standard” since 1965

ACTIVEAID, Inc. 101 Activeaid Road, Box 359 Redwood Falls, Minnesota 56283-0359, U.S.A.
Call: 507-644-2951 Toll Free: 1-800-533-5330 FAX: 507-644-2468 www.activeaid.com